InActiv Blue®

Pathogen inactivating and DNA/RNA stabilizing transport medium

document reference: FP09 I117 R01 B.3
update: 03/10/2024

INTENDED USE

InActiv Blue® is a pathogen inactivating and DNA/RNA stabilizing medium used for transport of naso- or oropharyngeal swabs, or saliva sampled from patients for in vitro diagnostic testing on the presence of genetic material of RNA/DNA by RT-qPCR or other molecular method.

MATERIAL INCLUDED

REF IB_TUB: PP tube (flat bottom) filled with 2 mL InActiv Blue®
A certificate of analysis and MSDS are available upon request or can be downloaded from the website www.inactivblue.com.

MATERIAL NEEDED (NOT PROVIDED)

- Nasopharyngeal or oropharyngeal swabs: it is advised to only use synthetic (flocked) fiber swabs with plastic or wire shafts and not to use calcium alginate or cotton swabs or swabs with wooden shafts as these may inhibit molecular tests¹. Ideally, the breakpoint of the swab is not higher than 8 cm (3.15 inch).
- Accessories for saliva collection (e.g. Blue Collect™).

SPECIMEN TYPE

- Nose/throat respiratory secretion collected by naso- or oropharyngeal swab.
- Saliva collected by spitting.

PERFORMANCE PARAMETERS

- InActiv Blue® completely inactivates pathogens within 5 minutes. More info upon request.
- InActiv Blue® maintains stability of RNA/DNA in saliva and swab for 30 days at 2-25 °C and 8 days at 37 °C.
- Proven stability after repeated freezing/thawing the sample.
- More details: www.inactivblue.com (download 'performance whitepaper' and 'Saliva study'²).

PATHOGEN INACTIVATION PERFORMANCE

- InActiv Blue® effectively inactivates a broad range of pathogens. Test results are shown in the table below.

Table: Inactivation status of pathogens when 60 minutes exposed to InActiv Blue®. [View "Table 1"](#)

STORAGE/DISPOSAL

- Store reagents between 2-25 °C.
- Keep away from direct (sun) light.
- Do not use after expiry date.
- InActiv Blue® needs to be disposed in accordance with local regulations for disposal of medical devices / hazardous substances.
- As indicated in the warning section below, contact between InActiv Blue® and bleach must be avoided.

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DO NOT DECONTAMINATE TUBES OR LAB ENVIRONMENT WITH BLEACH! InActiv Blue® contains guanidine thiocyanate that in combination with hypochlorite may form harmful cyanide gas.

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Harmful if swallowed, causes severe skin burns and eye damage, harmful to aquatic life with long lasting effects.

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All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting viruses. Always wear protective clothing when handling specimens and reagent (gloves, lab vest, surgical mask, eye/face protection).

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
- Make sure that you use the required personal protection. Follow the procedure for droplet protection: use gloves, coat, surgical mask and eye protection (goggles or face shield). Identify the sample.

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
- Collect the specimen with a suitable swab according to the instructions of the manufacturer of the swab. DO NOT PRE-WET the swab before touching the patient.

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
- While holding the swab, remove the cap from the test tube. Insert the swab into the test tube opening without touching the (outer) surface of the tube. Once the specimen is in the medium, spin around the swab.

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
- Break off the swab shaft (or alternatively cut the swab handle) so that the break point level corresponds with the test tube opening. Make sure that the breakpoint of the swab is not higher than the tube opening.

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- Hermetically close the tube by the screw cap. The sample is now stable for transport and/or storage prior to further processing in the laboratory.

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BIBLIOGRAPHY LITERATUR

- Centers for Disease Control and Prevention. (2020, February 14). Coronavirus Disease 2019 (COVID-19): Guidelines for Clinical Specimens. Retrieved February 27, 2020, from <https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>
- Jonckheere I, Faes L, Overmeire Y, De Vleeschauwer A, Vanden Daele L, Van Bruaene N, Vandecandelaere I, Merlaen B, van Cann J, Vandesompele J. 2021. 'Equivalence of saliva RT-qPCR testing to nasal-throat/nasopharyngeal swab testing in the general practitioner's setting to detect SARS-CoV-2'. medRxiv., Sept; 2021.09.30.21264181; Accepted for publication in J of Pediatrics, Perinatology and Child Health

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TECHNICAL SUPPORT TECHNISCHER SUPPORT

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