InActiv Blue®

Pathogen inactivating and DNA/RNA stabilizing transport medium

document reference: FP09 I117 R01 B.1
update: 20/09/2023

*For in vitro diagnostic use only
For professional use only*

INTENDED USE

InActiv Blue® is a pathogen inactivating and DNA/RNA stabilizing medium used for transport of naso- or oropharyngeal swabs, or saliva sampled from patients for in vitro diagnostic testing on the presence of genetic material of a virus by RT-qPCR or other molecular method (SARS-CoV-2, influenza A/B, RSV A/B).

GENERAL INFORMATION / COMPOSITION

InActiv Blue® abrogates the infectious potential of collected patient material without affecting the integrity of the RNA to be used in a diagnostic test based on RT-qPCR or other molecular method. InActiv Blue® also functions as a lysis buffer. It contains guanidine thiocyanate and N-lauroylsarcosine for pathogen and DNase/RNase inactivation. Complete product composition is provided in the MSDS

MATERIAL INCLUDED

REF IB_TUB: PP tube (flat bottom) filled with 2 mL InActiv Blue®
A certificate of analysis and MSDS are available upon request or can be downloaded from the website www.inactivblue.com.

MATERIAL NEEDED (NOT PROVIDED)

- Nasopharyngeal or oropharyngeal swabs: it is advised to only use synthetic (flocked) fiber swabs with plastic or wire shafts and not to use calcium alginate or cotton swabs or swabs with wooden shafts as these may inhibit molecular tests¹. Ideally, the breakpoint of the swab is not higher than 8 cm (3.15 inch).
- Accessories for saliva collection (e.g. Blue Collect™).

SPECIMEN TYPE

- Nose/throat respiratory secretion collected by naso- or oropharyngeal swab.
- Saliva collected by spitting.

PERFORMANCE PARAMETERS

- InActiv Blue® is able to completely inactivate SARS-CoV-2 virus in 1 minute.
- Also excellent inactivation of other viruses (e.g. RSV and influenza), yeast and bacteria. More info upon request.
- InActiv Blue® maintains the stability of SARS-CoV-2 RNA for at least 30 days (2-25 °C).
- Robust formulation: proven stability after repeated freezing/thawing the sample and 8 days exposure at 37 °C.
- More details: www.inactivblue.com (download 'performance whitepaper' and 'Saliva study'²).

STORAGE/DISPOSAL

- Store reagents/samples between 2 °C and 25 °C.
- Keep away from direct (sun) light.
- Do not use after expiry date.
- InActiv Blue® needs to be disposed in accordance with local regulations for disposal of medical devices / hazardous substances.
- As indicated in the warning section below, contact between InActiv Blue® and bleach must be avoided.

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DO NOT DECONTAMINATE TUBES OR LAB ENVIRONMENT WITH BLEACH! InActiv Blue® contains guanidine thiocyanate that in combination with hypochlorite may form harmful cyanide gas.

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Harmful if swallowed, causes severe skin burns and eye damage, harmful to aquatic life with long lasting effects.

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All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting viruses. Always wear protective clothing when handling specimens and reagent (gloves, lab vest, surgical mask, eye/face protection).

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Samples should be analyzed within 30 days upon sample collection. Invert tube 10 times to ensure that the complete inner surface of the tube (and swab) have been in contact with the inactivating buffer solution. Use validated laboratory procedures for PCR testing. Demonstrated compatibility with the test platforms listed in the table below (non-exhaustive list).

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instrument	RNA extraction	RT-qPCR
Abbott Alinity M	Alinity m Sample Prep Kit 1 #09N18-001	Alinity m SARS-CoV-2 AMP kit #09N78-090
CFX96 (Bio-Rad)	Real-Prep Viral DNA/RNA kit (BioSewoom)	STANDARD M nCoV Real-Time Detection kit (SD BIOSENSOR #11NCO10)
CFX96 Deep Well Real Time PCR detection System (Bio-Rad) Maelstrom 9600 (Tanbead)	Tanbead Nucleic Acid Extraction Kit (96) #W665A10	Allplex TM 2019 nCov Assay RP4520D59
CFX96 Deep Well Real Time PCR detection System (Bio-Rad) Microlab STARlet IVD (Seegene/Hamilton)	STARMag96x4 viral DNA/RNA 200c kit #EX00013C	Allplex TM 2019 nCov Assay RP4520D59
CFX384 (Bio-Rad #1855485)	Zymo Research's Quick-RNA 96 #R1053	Bio-Rad iTaq one-step RT-qPCR mix #1725141
	Norgen Biotek's Total RNA Purification 96-well Kit #24370	Bio-Rad iTaq one-step RT-qPCR mix #1725141
	magtivio's MagSi-NA Pathogens #MDKT00210960	Bio-Rad iTaq one-step RT-qPCR mix #1725141
Chemagic 360 Janus qPCR Janus reformatter (G3) Perkin Elmer QuantStudio7flex	Chemagic Viral DNA/RNA Kit special H96 #CMG-1033-S	Sars-CoV-2-RT-qPCR Reagent kit #3501-0010
cobas Liat PCR system (Roche)	COBAS SARS-CoV-2	COBAS SARS-CoV-2
cobas 6800 system (Roche)	cobas omni reagent for sample preparation	cobas SARS-CoV-2: #09175431190 cobas Utility Channel with probes and primers for SARS -CoV-2 E gene
EMAG (Biomérieux) Stratagene MX3000-3005 (Agilent)	NUCLISENS easyMAG (Biomérieux #280130-35; #280146; #200292)	Superscript III Platinum One-Step qRT-PCR system (Invitrogen #11732-088)
Hamilton Starlet	96X4 viral DNA/RNA 200 C kit #EX00013C	Allplex SARS-CoV-2 #RV10248X
Hamilton Starlet/Tanbead	OptiPure Viral Auto Plate #W665A10	Allplex SARS-CoV-2 #RV10248X
Lightcycler 480 (Roche)	MagnaPure 96 and Viral NA Small volume kit 0654388001	LightMix Modular Sarbecovirus SARS-CoV2 PCR kit #50-0776-96
Kingfisher FLEX (Thermofisher) Quantstudio 5 (Thermofisher)	Thermofisher's MagMAX Viral/ Pathogen II kit #A48383	TaqPath COVID-19 CE-IVD RT-PCR kit #A48067
Kingfisher FLEX (Thermofisher) Quantstudio 7 (Thermofisher)	Thermofisher's MagMAX Viral/ Pathogen II kit #A48383	TaqPath COVID-19 CE-IVD RT-PCR kit #A48067
Lightcycler 96 (Roche)	magtivio's MagSi-NA Pathogens #MDKT00210960	Bio-Rad iTaq one-step RT-qPCR mix #1725141
BioFire FilmArray (Biomérieux)	BioFire Respiratory 2.1 plus Panel (Biomérieux)	

BIBLIOGRAPHY LITERATUR

- Centers for Disease Control and Prevention. (2020, February 14). Coronavirus Disease 2019 (COVID-19): Guidelines for Clinical Specimens. Retrieved February 27, 2020, from https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html
- Jonckheere I, Faes L, Overmeire Y, De Vleeschauer A, Vanden Daele L, Van Bruaene N, Vandecandelaere I, Merlaen B, van Cann J, Vandesompele J. 2021. 'Equivalence of saliva RT-qPCR testing to nasal-throat/nasopharyngeal swab testing in the general practitioner's setting to detect SARS-CoV-2'. medRxiv., Sept; 2021.09.30.21264181; Accepted for publication in J of Pediatrics, Perinatology and Child Health

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
1 Make sure that you use the required personal protection. Follow the procedure for droplet protection: use gloves, coat, surgical mask and eye protection (goggles or face shield). Identify the sample.

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
2 Collect the specimen with a suitable swab according to the instructions of the manufacturer of the swab. DO NOT PRE-WET the swab before touching the patient.

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
3 While holding the swab, remove the cap from the test tube. Insert the swab into the test tube opening without touching the (outer) surface of the tube. Once the specimen is in the medium, spin around the swab.

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
4 Break off the swab shaft (or alternatively cut the swab handle) so that the break point level corresponds with the test tube opening. Make sure that the breakpoint of the swab is not higher than the tube opening.

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5 Hermetically close the tube by the screw cap. The sample is now stable for transport and/or storage prior to further processing in the laboratory.

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