










## InActiv Blue® - Best in Class

*InActiv Blue® is a CE marked, virus inactivating and RNA stabilizing transport medium for molecular analyses including RT-qPCR*



-  Perfect for preservation of nose/throat swabs and saliva
-  Complete inactivation of SARS-CoV-2 **after 5 minutes**
-  Safe and stable preservation of RNA up to 30 days
-  No cold chain required
-  CE marked
-  Blue color as visual pipetting control
-  Independent validation
-  Production and testing according to ISO 13485 certified labs
-  Developed in Belgium

### Products

### Order code

Virus transport medium and swab (CE marked)

|  |               |
|--|---------------|
| InActiv Blue® tube (1x PET tube (100 x 13 mm) with 2 ml virus inactivating and RNA stabilizing transport medium) | <b>456604</b> |
|--|---------------|

|  |                |
|--|----------------|
| Swab, nasopharyngeal, flocked nylon tip, sterile, length 151 mm, breakpoint at 80 mm | <b>4800357</b> |
|--|----------------|

|                       |               |
|-----------------------|---------------|
| Saliva Collection Kit | <b>IB_COL</b> |
|-----------------------|---------------|

For RESEARCH ONLY

|   |              |
|---|--------------|
| InActiv Blue® (1x PETG bottle - 250 ml) | <b>IB250</b> |
|---|--------------|

|  |               |
|--|---------------|
| InActiv Blue® (1x PETG bottle - 1000 ml) | <b>IB1000</b> |
|--|---------------|

### LEARN MORE ABOUT INACTIV BLUE?

Mail to: [info@inactivblue.com](mailto:info@inactivblue.com)

Website: [www.inactivblue.com](http://www.inactivblue.com)

We are happy to answer all your questions!